

NEX IODIO P2

REFS & PACKAGING



REFS

code: **NEXP2EN001B:**
ENGLISH LABEL, 6X40 UNITS, WITH DISPENSER BOX 40 UNITS, WITH NAILCLEANER

code: **NEXP2EN005B:**
ENGLISH LABEL, 240 UNITS, BULK, WITH NAILCLEANER

code: **NEXP2EN009B:**
ENGLISH LABEL, 240 UNITS, BULK, WITHOUT NAILCLEANER

code: **NEXP2EN013B:**
ENGLISH LABEL, 6X40 UNITS, WITH DISPENSER BOX 40 UNITS, WITHOUT NAILCLEANER

code: **NEXP2EN017B:**
ENGLISH LABEL, 100 UNITS, CARTON BOX WITH SIDE OPENING, WITH NAILCLEANER

code: **NEXP2EN021B:**
ENGLISH LABEL, 100 UNITS, CARTON BOX WITH SIDE OPENING, WITHOUT NAILCLEANER

code: **NEXP2SP025B:**
SPANISH LABEL, 6X40 UNITS, WITH DISPENSER BOX 40 UNITS, WITH NAILCLEANER

code: **NEXP2SP027B:**
SPANISH LABEL, 240 UNITS, BULK, WITH NAILCLEANER

code: **NEXP2SP029B:**
SPANISH LABEL, 240 UNITS, BULK, WITHOUT NAILCLEANER

code: **NEXP2SP031B:**
SPANISH LABEL, 6X40 UNITS, WITH DISPENSER BOX 40 UNITS, WITHOUT NAILCLEANER

code: **NEXP2SP033B:**
SPANISH LABEL, 100 UNITS, CARTON BOX WITH SIDE OPENING, WITH NAILCLEANER

code: **NEXP2SP035B:**
SPANISH LABEL, 100 UNITS, CARTON BOX WITH SIDE OPENING, WITH NAILCLEANER

PACKAGING

DISPENSER BOX 40 UNITS:
145mmx130mmx540mm (H)

CARTON BOX 240 UNITS:
40cmx31cmx56cm (H)

CARTON BOX 100 UNITS:
40cmx31cmx27cm (H)

PALLET CARTON BOX 240 UNITS:
5.760units (H 1,85m) / 7.680 units (H 2,30m)

PALLET CARTON BOX 100 UNITS:
4.800units (H 1,85m) / 6.400 units (H 2,30m)


NEX
MEDICAL Surgical scrub

NEX IODIO P2[®] SURGICAL SCRUB

Hand Hygiene is the basis of a correct and effective infection prevention in the hospital environment.

Surgical hand antiseptics, to destroy transient micro-organisms and inhibit the growth of resident micro-organisms, is routinely carried out before undertaking invasive procedures.

Antisepsis may reduce the risk of surgical site infections in patients.



NEX IODIO P2 PVP-I 7,5% surgical brush/sponge

Description

NEX IODIO P2 is a single use surgical hand scrub/brush with nailcleanser impregnated with Povidone Iodine antiseptic solution.

Indications

Antiseptic surgical scrub. Disinfection of the hands prior to surgery procedures or in the operating theatre.

Formulation Composition

Ingredients for 100g.:
Povidone iodine g. 7,5
Cofomulants and water q.s. to g. 100
free iodine 0,75%

Microbiological & skin irritation tests

Bactericidal:
EN 1040, pr-EN 12054 (Surgical handrub and handwash) and pr-EN 12791 (Surgical hand disinfection)
Dermal irritation test and sensitization
(OECD guidelines 404, 410, ISO 10993)

Materials

Brush: PE medical grade
Sponge: Polyurethane
Nailcleanser: PP medical grade
All components are LATEX free

Manufacturer

NEX MEDICAL
20010 Casorezzo (MI) - ITALY
reg. PMC N. 26/2011/Off.299 / UNI EN ISO 9001:2008 / UNI EN ISO 13485:2004

Product Ref: NEX IODIO P2

Colour code: brown



Ergonomic shape to allow a perfect grip
High flexibility to ensure accurate cleaning
Soft and kind to the skin
With nailcleanser for accurate cleaning under the nails
Advanced bristles system
Easy peel-off single wrapping
Single use

NEX MEDICAL Surgical scrub

www.nexmedical.com

sales@nexmedical.it

**NEX IODIO P2
SURGICAL SCRUB**

Surgical brush/sponge with nailcleanser - Povidone iodine 7,5% (0,75% available iodine)
NEX IODIO P2 is a single use surgical brush/sponge impregnated with approximately 20ml of Povidone iodine antiseptic solution

Uses: surgical hand scrub, for cleansing hands and arms prior to surgery
Active ingredients: Povidone iodine 7,5% (equivalent to 0,75% available iodine)
Direct action towards bacteria: gram+, gram- EN 12791, EN 12054, EN 1040
Directions: wet hands and forearms with water. Scrub for at least 3 minutes paying particular attention to the nails, cuticles and interdigital spaces. Rinse and repeat scrub. Discard brush-sponge. Rinse hands and forearms thoroughly and dry with sterile cloth.

Warnings: For external use only. (R36) Irritating to eyes (S2) Keep out of reach of children (S28) Avoid contact with eyes and ears (S28) In case of contact with eyes, rinse immediately with plenty of water and seek medical advice (S46) If swallowed, seek medical advice immediately and show this container or label

Ingredients 100 g contain: Povidone-iodine g. 7,5 ; Cofomulants and water q.s.to g. 100

Italian Ministry of Health registration PMC n°18975 (D.G.F.D.M. 1405) Store at room temperature
Manufacturer's Address: Nex Medical s.a.s. - 20010 Casorezzo (MI) - ITALY Latex free
Made in Italy www.nexmedical.it info@nexmedical.it

NEX CLOREX C2

REFS & PACKAGING



REFS

code: **NEXC2EN002B:**
ENGLISH LABEL, 6X40 UNITS, WITH DISPENSER BOX 40 UNITS, WITH NAILCLEANER

code: **NEXC2EN006B:**
ENGLISH LABEL, 240 UNITS, BULK, WITH NAILCLEANER

code: **NEXC2EN010B:**
ENGLISH LABEL, 240 UNITS, BULK, WITHOUT NAILCLEANER

code: **NEXC2EN014B:**
ENGLISH LABEL, 6X40 UNITS, WITH DISPENSER BOX 40 UNITS, WITHOUT NAILCLEANER

code: **NEXC2EN018B:**
ENGLISH LABEL, 100 UNITS, CARTON BOX WITH SIDE OPENING, WITH NAILCLEANER

code: **NEXC2EN022B:**
ENGLISH LABEL, 100 UNITS, CARTON BOX WITH SIDE OPENING, WITHOUT NAILCLEANER

code: **NEXC2SP026B:**
SPANISH LABEL, 6X40 UNITS, WITH DISPENSER BOX 40 UNITS, WITH NAILCLEANER

code: **NEXC2SP028B:**
SPANISH LABEL, 240 UNITS, BULK, WITH NAILCLEANER

code: **NEXC2SP030B:**
SPANISH LABEL, 240 UNITS, BULK, WITHOUT NAILCLEANER

code: **NEXC2SP032B:**
SPANISH LABEL, 6X40 UNITS, WITH DISPENSER BOX 40 UNITS, WITHOUT NAILCLEANER

code: **NEXC2SP034B:**
SPANISH LABEL, 100 UNITS, CARTON BOX WITH SIDE OPENING, WITH NAILCLEANER

code: **NEXC2SP036B:**
SPANISH LABEL, 100 UNITS, CARTON BOX WITH SIDE OPENING, WITH NAILCLEANER

PACKAGING

DISPENSER BOX 40 UNITS:
145mmx130mmx540mm (H)

CARTON BOX 240 UNITS:
40cmx31cmx56cm (H)

CARTON BOX 100 UNITS:
40cmx31cmx27cm (H)

PALLET CARTON BOX 240 UNITS:
5.760units (H 1,85m) / 7.680 units (H 2,30m)

PALLET CARTON BOX 100 UNITS:
4.800units (H 1,85m) / 6.400 units (H 2,30m)

NEX Surgical scrub
MEDICAL

NEX CLOREX C2[®] SURGICAL[®] SCRUB

Hand Hygiene is the basis of a correct and effective infection prevention in the hospital environment.

Surgical hand antiseptics, to destroy transient micro-organisms and inhibit the growth of resident micro-organisms, is routinely carried out before undertaking invasive procedures.

Antisepsis may reduce the risk of surgical site infections in patients.



NEX CLOREX C2 CHG 4% surgical brush/sponge

Description

NEX CLOREX C2 is a single use surgical hand scrub/brush with nailcleaner impregnated with Chlorhexidine gluconate antiseptic solution.

Indications

Antiseptic surgical scrub. Disinfection of the hands prior to surgery procedures or in the operating theatre.

Formulation Composition

Ingredients for 100g.:
Chlorhexidine gluconate g. 4,0
Coformulants and water q.s. to g. 100

Microbiological & skin irritation tests

Bactericidal:
EN 1040, pr-EN 12054 (Surgical handrub and handwash) and pr-EN 12791 (Surgical hand disinfection)
Dermal irritation test and sensitization
(OECD guidelines 404, 410, ISO 10993)


Materials

Brush: PE medical grade Sponge: Polyurethane
Nailcleaner: PP medical grade
All components are LATEX free

Manufacturer

NEX MEDICAL
20010 Casorezzo (MI) - ITALY
reg. PMC N. 26/2011/Off.299 / UNI EN ISO 9001:2008 / UNI EN ISO 13485:2004

Product Ref: NEX CLOREX C2
Colour code: orange



Ergonomic shape to allow a perfect grip
High flexibility to ensure accurate cleaning
Soft and kind to the skin
With nailcleaner for accurate cleaning under the nails
Advanced bristles system
Easy peel-off single wrapping
Single use

NEX
MEDICAL Surgical scrub

www.nexmedical.com

sales@nexmedical.it

NEX CLOREX C2[®] SURGICAL[®] SCRUB

Surgical brush/sponge with nailcleaner - Chlorhexidine digluconate 4%
NEX CLOREX C2 is a single use surgical brush/sponge impregnated with approximately 20ml of Chlorhexidine Gluconate antiseptic solution

Uses: surgical hand scrub, for cleansing hands and arms prior to surgery

Active ingredients: Chlorhexidine digluconate 4% (Antiseptic)

Direct action towards bacteria gram+, gram- EN 12791, EN 12054, EN 1040

Directions: wet hands and forearms with water. Scrub for at least 3 minutes paying particular attention to the nails, cuticles and interdigital spaces. Rinse and repeat scrub. Discard brush-sponge. Rinse hands and forearms thoroughly and dry with sterile cloth.

Warnings For external use only. Incompatible with Hydrogen Peroxide
(R36) Irritating to eyes (S2) Keep out of reach of children. (S25) Avoid contact with eyes and ears. (S26) In case of contact with eyes, rinse immediately with plenty of water and seek medical advice (S46) If swallowed, seek medical advice immediately and show this container or label

Ingredients 100 g contain : Chlorhexidine digluconate g. 4,0 ; Coformulants and water q.s.to g. 100

Italian Ministry of Health registration PMC n°18976 (D.G.F.D.M. 1400)

Manufacturer's Address: Nex Medical s.a.s. - 20019 Settimo Milanese (MI) - ITALY

Made in Italy

Reg. PMC N. 20/2008/Off.280

Store at room temperature

Latex free

www.nexmedical.it info@nexmedical.it

NEX D1 DRY

REFS & PACKAGING



REFS

code: **NEXD1EN004B:**
6X40 UNITS, WITH DISPENSER BOX 40 UNITS, WITH NAILCLEANER

code: **NEXD1EN016B:**
6X40 UNITS, WITH DISPENSER BOX 40 UNITS, WITHOUT NAILCLEANER

code: **NEXD1EN008B:**
240 UNITS, BULK, WITH NAILCLEANER

code: **NEXD1EN012B:**
240 UNITS, BULK, WITHOUT NAILCLEANER

code: **NEXD1EN020B:**
100 UNITS, CARTON BOX WITH SIDE OPENING, WITH NAILCLEANER

code: **NEXD1EN024B:**
100 UNITS, CARTON BOX WITH SIDE OPENING, WITHOUT NAILCLEANER

PACKAGING

DISPENSER BOX 40 UNITS:
145mmx130mmx540mm (H)

CARTON BOX 240 UNITS:
40cmx31cmx56cm (H)

CARTON BOX 100 UNITS:
40cmx31cmx27cm (H)

PALLET CARTON BOX 240 UNITS:
5.760units (H 1,85m) / 7.680 units (H 2,30m)

PALLET CARTON BOX 100 UNITS:
4.800units (H 1,85m) / 6.400 units (H 2,30m)

NEX Surgical scrub
MEDICAL

NEX D1 DRY[®] SURGICAL SCRUB

Hand Hygiene is the basis of a correct and effective infection prevention in the hospital environment.

Surgical hand antisepsis, to destroy transient micro-organisms and inhibit the growth of resident micro-organisms, is routinely carried out before undertaking invasive procedures.

Antisepsis may reduce the risk of surgical site infections in patients.



NEX D1 DRY - sterile EO surgical brush/sponge

Description

NEX D1 DRY is a single use brush/sponge with nailcleaner
Without DETERGENT
Sterilized by ETHYLENE OXIDE (EO)
Single use
Expiration: 5 years

Indications

Surgical scrubbing of hands for antiseptic pre-op procedures.
NEX D1 DRY is compatible with any detergent or antiseptic cleansing solution.

Materials

Brush: PE medical grade
Sponge: polyurethane
Nailcleaner: PP medical grade
All components are LATEX free

Microbiological & skin irritation tests

Hemolysis test (ISO 10993-4:2008)
Cytotoxicity test (ISO 10993-5:2009)
Skin irritation test (ISO 10993-10:2002 / Amd 1:2006)

Manufacturer

NEX MEDICAL
20010 Casorezzo (MI) - ITALY
reg. PMC N. 26/2011/Off.299 / UNI EN ISO 9001:2008 / UNI EN ISO 13485:2004

Product Ref: NEX D1 DRY
Colour code: green




Ergonomic shape to allow a perfect grip
High flexibility to ensure accurate cleaning
Soft and kind to the skin
With nailcleaner for accurate cleaning under the nails
Advanced bristles system
Easy peel-off single wrapping
Single use

NEX
MEDICAL Surgical scrub

www.nexmedical.com

sales@nexmedical.it



NEX SPONGE LINE

ready to use surgical pre-op & prep sponges

Indications

NEX SPONGE is a complete line of single use surgical sponges impregnated with Chlorhexidine digluconate or PVP-I antiseptic solutions, ready to use.

Description

Pre-operative skin preparation.

Materials

Sponge: Polyurethane
All components are LATEX free

Pre-operative skin preparation & disinfection (aqueous solutions)



NEX CLOREX C2 CHG 2% SPONGE

ACTIVE: Chlorhexidine digluconate 2%

Composition 100g: Chlorhexidine gluconate 2,0 g, coformulants and water up to 100g

Available in two sizes: small (9cmx4cmx2,5cm) and medium (12cmx8cmx2,5cm)



NEX CLOREX C2 CHG 4% SPONGE

ACTIVE: Chlorhexidine digluconate 4%

Composition 100g: Chlorhexidine gluconate 4,0 g, coformulants and water up to 100g

Available in two sizes: small (9cmx4cmx2,5cm) and medium (12cmx8cmx2,5cm)



NEX CLOREX P2 PVP-I 1% SPONGE

ACTIVE: Povidone Iodine at 10% (1% free iodine)

Composition 100g: Povidone Iodine at 10% 10,0 g, coformulants and water up to 100g

Available in two sizes: small (9cmx4cmx2,5cm) and medium (12cmx8cmx2,5cm)

Pre-operative disinfection & washing (soapy solutions)



NEX CLOREX C2 CHG 2% WASHING SPONGE

ACTIVE: Chlorhexidine digluconate 4%

Composition 100g: Chlorhexidine gluconate 4,0 g, coformulants and water up to 100g

Size: 12cmx8cmx2,5cm



NEX CLOREX C2 CHG 4% WASHING SPONGE

ACTIVE: Chlorhexidine digluconate 4%

Composition 100g: Chlorhexidine gluconate 4,0 g, coformulants and water up to 100g

Size: 12cmx8cmx2,5cm

NEX
MEDICAL Surgical scrub